

## Background

The prognosis and treatment of patients with resected colon cancer are based on the TNM staging classification. However, this approach fails to account for the important role of the immune system.

The immune response can now be accurately measured in clinical practice with the Immunoscore® Colon assay. This assay provides a standardized quantification of CD3+ and CD8+ T cell densities both in the tumor center and invasive margin.


Immunoscore® was shown to predict patient outcome more accurately than TNM classification plus other parameters including microsatellite instability status in a large international multicenter study, supporting its implementation as a new component of a TNM-immune classification of colon cancer.

The objective of PROSCORE study is to assess the impact of Immunoscore® on adjuvant chemotherapy decision making in patients with resected stage II-III colon cancer.

## Study design

PROSCORE is a multicenter, non-randomized, prospective study.

The medical decisions (MDs) will be taken during one multidisciplinary meeting (MM). Alternatively, the MDs can be taken during two distinct MMs.



## Objectives

### Primary objective :

To assess whether the Immunoscore® (IS) Colon test results will impact the choice of therapy, during a multidisciplinary meeting, in two cohorts of patients (stage II and stage III colon cancer) after curative-intent surgery of non-metastatic colon cancer for which adjuvant chemotherapy is being considered.

### Secondary objectives :

To assess:

- The consistency of decisions with the current recommendations
- The decision-making impact of IS Colon test results on treatment selection according to clinical sub-groups
- The feasibility of IS Colon test in clinical standard practice, samples and information flow
- The correlation of IS Colon test result with relapse-free survival (RFS), disease-free survival (DFS) and disease-specific survival (DSS)
- RFS rate at 3 and 5 years
- Overall survival at 5 years
- The budget impact of IS Colon test use

## Trial Rationale

TNM classification in combination with conventional clinicopathological factors provides the most reliable reference for routine prognostication and for guiding treatment decisions for colon carcinoma. However, the information provided by those classification tools remains imperfect in predicting the outcome of patients

Immunoscore® allows us to identify sub-groups of patients (see fig.1 &2):

- 1) **Patients at high risk of relapse** who may benefit from adjuvant chemotherapy (for stage II patients) or treatment intensification (for stage III patients).
- 2) **Patients at low risk of relapse** who may benefit from reduced adjuvant chemotherapy (decreased intensity and/or shorter duration) in stage III patients or absence of adjuvant chemotherapy in stage II patients.

Immunoscore® is an additional objective criterion for making individualized therapeutic decision with an improved risk assessment. The use of Immunoscore® aims to improve the treatments for patients and potential reduce treatment costs for the national health system.

## IS test result disclosing

During a multidisciplinary meeting, a **first therapeutic decision (n°1)** concerning adjuvant chemotherapy will be made and documented.

Once the initial decision is made, the IS Colon test result will be **disclosed** (eCRF and IVRS systems).

A second **therapeutic decision (n°2)** will then be made and documented.

## Statistical considerations

The primary endpoint is the **rate of therapeutic modifications** due to IS Colon test result.

Sample size: **280 patients enrolled** (140 stage II and 140 stage III CRC)

With a **one-stage Fleming design** ( $\alpha=5\%$ ,  $\beta=5\%$ ,  $p_0=10\%$  and  $p_1=20\%$ ) we will include 133 evaluable per cohort. Taking into account an estimate of 5% of non-evaluable patients, we will enrol a total of 140 patients in each cohort.

The strategy will be considered effective if we observed at least **19 successes out of 133 evaluable patients**. A success being defined as a therapeutic modification due to the IS Colon test result.




Figure 1 - TTR in stage II colon cancer



Classification	Events/Total	2 Year KM Est.(95% CI)	Hazard Ratio (95% CI)
Immunoscore High (n=134)	23/134	91.8 (87.3-96.6%)	0.55 (0.32-0.92)
Immunoscore Low (n=153)	43/153	77.5 (71.3-84.4%)	Reference

Figure 2 - DFS in stage III colon cancer

### References:

Pages, F et al., The Lancet, May 2018  
Sinicrope, F et al., J Clin Oncol 36, 2018

**Study start: January 2019**  
**End of recruitment: June 2020**

### Conflict of interests

First author disclosed relationships with Amgen, Roche, Sanofi-Aventis, Servier, Merck, Merck Serono, HaloDx, Bayer and Shire.  
Co-authors F. Hermitte and A. Catteau disclosed relationships with HaloDx

### Information

Please contact Dr. David MALKA, study coordinator: david.malka@gustaveroussy.fr